

Freestyle Edge Firmware Release Notes – Document Version July 10, 2020

Version 1.0.536. Released for beta testing on 7/10/20.

Bug Fixed: Macro playback issues created in v.1.0.535

Version 1.0.535. Released for beta testing on 5/1/20.

New Feature: Ability to navigate layers with a held modifier key

Bug Fixed: Media actions in macros

Version 1.0.532. Released for beta testing on 8/5/19.

Bug Fixed: Tap-and-Hold performance improved

Version 1.0.510. Released for beta testing on 5/29/19.

Bug Fixed: Tap-and-Hold performance improved

Version 1.0.480. Released for beta testing on 5/3/19. (Download the latest version of the SmartSet App)

New Feature: Tap-and-Hold

- Use the token [t&hxxx] to assign two unique actions to a single key based on how long the key is depressed (where xxx is the “timing delay” from 1ms to 999ms). The Tap action occurs when the key is pressed and released faster than the timing delay, and the Hold action occurs when the key is pressed and held longer than the timing delay.
 - Examples:
 - Caps Lock becomes Control when held for longer than 250ms
 - Syntax: [caps]>[caps][t&h250][lctrl]
 - Left Shift becomes Escape when tapped & released in less than 100ms
 - Syntax: [lshift]>[esc][t&h100][lshift]
 - *Notes
 - Not suitable for assigning dual actions to primary typing keys or keys that have a shifted action
 - Not all key actions supported

New Feature: Multimodifiers “Hyper” and “Meh”

- Use the action tokens [hyper] and [meh] in a layout file to send a multi-modifier command (from any key) that can be combined with another key for use with custom shortcuts via 3rd party software.
 - Hyper = Left Shift + Left Control + Left Alt + Left Windows
 - Syntax: [position token]>[hyper]
 - Meh = Left Shift + Left Control + Left Alt
 - Syntax: [position token]>[meh]

Version 1.0.372. Released for beta testing on 4/29/19.

Bug Fixed: Breathe mode re-enabled

Version 1.0.340. Released for beta testing on 3/13/18.

New Feature: Macro capacity increase to 100 macros (7,200 total characters)

Version 1.0.101. Released for beta testing on 1/11/19.

Bug Fixed: Game Mode no longer blocking F4

Bug Fixed: Embedded keys can trigger macros

Version 1.0.90. Released for beta testing on 12/16/17.

New Feature: Random & Custom Macro delays added (now supported in SmartSet App)

Version 1.0.74. Released for beta testing on 11/8/17.

New Feature: Added HID Code Support (*Direct Programming required)

New Feature: FAT12 Reformat shortcut (SmartSet + Shift + Ctrl + F)

Version 1.0.64. Released for beta testing on 10/13/17.

Bug Fixed: Improved backlighting brightness range

Bug Fixed: Token for fn [kp0] fixed